



Patient Name: _____

Reading ID: _____

This Oximetry test has been ordered by your physician to monitor your blood Oxygen levels while asleep. In order to obtain a quality test, it is essential that you closely follow these instructions. If you have any questions, you can contact VirtuOx at (877) 337-7111 for further assistance. **You must stay awake for the first 10 minutes of the test!**

VirtuOx is the testing facility that will process the results of the test you will be taking and forward them to your physician. Per Medicare guidelines, the company that delivers the testing device to you cannot assist with the performance or instruction of the test and can only act as the courier.

The DME that couriered the device is: _____ **Phone:** _____

Your physician would like you to test under the following **Test Conditions:**

- | | |
|--|---|
| <input type="checkbox"/> Overnight on Room Air | <input type="checkbox"/> Overnight on APAP / CPAP / BiPAP |
| <input type="checkbox"/> Overnight on Oxygen @ _____ LPM | <input type="checkbox"/> Overnight on Dental Device |

Step 1: Follow Test Condition above

It is suggested that you remove fingernail polish and coverings from the finger you will be using.

Step 2: Sensor should be attached as shown below to device



Step 3: Attach sensor to Index Finger, apply tape or Band Aid to secure sensor to finger



Step 4: Turn on Device
Numbers will display



Step 5: Go to sleep, wear all night



Step 6: When you wake up, turn off device and call DME above for pick up



IF THE SENSOR FALLS OFF DURING THE NIGHT:

- 1) Replace the sensor on your finger
- 2) Check that the oximeter is on and displaying numbers.
- 3) If the oximeter is OFF, press the power button to turn it on.

IF THE INDICATOR LIGHT BLINKS YELLOW OR RED FOR SEVERAL MINUTES: 1) Reposition the sensor on finger. 2) Continue the test.

IF THE LOW BATTERY INDICATOR IS YELLOW: 1) **DO NOT REPLACE BATTERIES.** 2) Continue the test. 3) Advise the DME when they pick up the equipment.