## **CHEYNE-STOKES RESPIRATION**

# Validation of ApneaLink as Screening Device for Cheyne-Stokes Respiration

Gerhard Weinreich, PhD<sup>1</sup>; Jeff Armitstead, PhD<sup>2</sup>; Volker Töpfer, MD<sup>1</sup>; You-Ming Wang, MD<sup>1</sup>; Yi Wang, MD<sup>1</sup>; Helmut Teschler, MD<sup>1,2</sup>

<sup>1</sup>Ruhrlandklinik, Department of Respiratory Medicine, Faculty of Medicine, University of Duisburg-Essen, Essen, Germany; <sup>2</sup>University of New South Wales, Graduate School of Biomedical Engineering, Sydney, Australia

**Study Objectives:** Polysomnography (PSG) is the gold standard in diagnosing sleep disordered breathing (SDB). Diagnostic analysis of SDB is time-consuming and labor-intensive and entails long waiting lists for patients. The aim of this study was to assess the validity of a flow-based screening classifier as an automated diagnostic test for Cheyne-Stokes respiration (CSR).

Setting: Sleep laboratory

**Participants:** 70 study subjects (28 with obstructive sleep apnea [OSA], 20 with CSR, 11 with CSR+OSA and 11 without SDB)

**Measurements:** The nasal cannula flow signal was analyzed by ApneaLink (ResMed, Sydney, Australia), based on a classifier algorithm using pattern recognition. In a simultaneous PSG, results were compared with manual scoring of respiratory events by certified sleep experts. **Results:** For detecting CSR we obtained a sensitivity of 87.1% (95%)

CHEYNE-STOKES RESPIRATION (CSR) IS AN ABNOR-MAL PERIODIC PATTERN OF BREATHING PREDOMI-NANTLY IN NREM SLEEP. CSR IS CHARACTERIZED by waxing and waning flow amplitude followed by apnea or hypopnea (Figure 1). This pattern is caused by instability in the feedback control involved in the chemical regulation of breathing. CSR affects 30% to 40% of patients with congestive heart failure<sup>1-4</sup> and can also be seen in patients with stroke.<sup>5</sup> CSR is associated with sleep disruption and decreased survival rates compared to subjects with a similar ejection fraction but without CSR.<sup>6-8</sup>

Diagnostic assessment of SDB is based on polysomnography (PSG) recorded during the night. Because diagnostic analysis is labor-intensive and time-consuming, patients must wait 2-10 months for a PSG.<sup>9</sup> To prioritize patients for PSG, it is desirable to have a screening device that can distinguish between the various manifestations of SDB. ApneaLink is a single-signal screening device that measures flow via a nasal cannula connected to a pressure transducer. Recent studies showed an excellent correlation of ApneaLink versus PSG regarding the apnea-hypopnea index (AHI).<sup>10,11</sup>

The newest version of ApneaLink provides a classifier for the auto-detection of CSR. The classifier of ApneaLink is based on pattern recognition, a concept which is also used for voice recognition and radar detection.<sup>12-15</sup> The nasal cannula flow signal of ApneaLink is used for pattern recognition and hence classification of CSR.

Submitted for publication September, 2008 Submitted in final revised form January, 2009 Accepted for publication January, 2009

Address correspondence to: Gerhard Weinreich, PhD, Ruhrlandklinik, Department of Respiratory Medicine, Faculty of Medicine, University of Duisburg-Essen, Essen, Germany; E-mail: gerhard.weinreich@uni-due.de

confidence interval 75.3% to 98.9%), a specificity of 94.9% (95% confidence interval 87.9% to 100%), a positive likelihood ratio of 17.0, and a negative likelihood ratio of 0.14. The area under the curve (AUC) of the according receiver operating characteristic (ROC) curve was 93.4%. This resulted in an accuracy of 91.4% for identifying CSR. **Conclusion:** In this study we demonstrated that the screening classifier was able to detect CSR with high diagnostic accuracy. Hence, ApneaLink equipped with CSR classifier is an appropriate screening tool which may help to prioritize patients with CSR for PSG. **Keywords:** ApneaLink, Cheyne-Stokes respiration, obstructive sleep apnea, pattern recognition, sleep disordered breathing, screening **Citation:** Weinreich G; Armitstead J; Töpfer V; Wang YM; Teschler H. Validation of apnealink as screening device for cheyne-stokes respiration. *SLEEP* 2009;32(4):553-557.

In recent years attempts have been made to automatically detect SDB. El-Solh et al<sup>16</sup> evaluated a neural network based on spectral analysis of arterial oxygen saturation SpO<sub>2</sub> in order to determine CSR. Our group analyzed a neural network relying on spectral entropy S of the airflow for identifying OSA and CSR.<sup>17</sup> The aim of the present study was to evaluate the diagnostic accuracy of the ApneaLink classifier for CSR against simultaneous PSG.

## METHODS

### Subjects

The study population consisted of 70 subjects referred to our sleep laboratory for PSG testing (Table 1). Firstly, we included 41 cardiology patients with congestive heart failure and high pre-test probability of CSR. Secondly, to test the device in a mixed population (CSR, OSA, CSR+OSA, no SDB), we included 29 subjects with high pre-test probability of OSA referred from general practitioners. An institutional review board approved the study, and subjects provided informed consent prior to participation in this study.

#### Polysomnography


Polysomnography was performed using a digital polygraph (Embla, Broomfield, USA) including 2 electroencephalographic (EEG) channels, 2 electrooculograms (EOG), submental electromyogram (EMG), tibialis EMG, rib cage and abdominal inductance pneumograms, pulse oximeter (Nonin 8500, Minnesota, USA), body position, and nasal cannula connected via a T-piece allowing simultaneous recording of flow signal by ApneaLink and the PSG system. Sleep staging and arousals were scored using 30-s epochs according to Rechtschaffen and 
 Table 1—Demographic Characteristics, SDB Classification, and Comorbidities for All 70 Subjects


Demographic Characteristic	Results		
Age	$60.3 \pm 11.5$ years		
Sex, f/m	15/55		
BMI	$27.5 \pm 4.4 \text{ kg/m}^2$		
OSA	28		
CSR	20		
OSA + CSR	11		
no SDB	11		
Congestive heart failure	41		
Hypertension	36		
Diabetes	13		
Asthma	9		
Chronic obstructive bronchitis	14		
Allergies	5		


Kales<sup>19</sup> and American Sleep Disorders Association.<sup>20</sup> An apnea was defined as cessation of flow for  $\geq 10$  s, and a hypopnea as reduction of  $\geq$  50% in flow amplitude lasting  $\geq$  10 s. The respiratory event was scored as obstructive apnea if it met apnea criteria and was associated with continued or increased respiratory effort throughout the entire period of absent airflow.<sup>21</sup> The respiratory event was scored as central apnea if it met apnea criteria and was associated with absent respiratory effort throughout the entire period of absent airflow. According to the definitions of the American Academy of Sleep Medicine<sup>21</sup> CSR was defined if i)  $\geq$  5 central apneas or hypopneas/h occurred, or ii) the cyclic crescendo and decrescendo change in breathing amplitude had a duration of  $\geq 10$  consecutive minutes. In PSG an apnea-hypopnea index (AHI) was calculated by dividing the total number of apneas and hypopneas by the total hours of sleep time recorded over the night.

#### ApneaLink

ApneaLink (ResMed, Sydney, Australia) is a single-signal screening device which measures flow via a nasal cannula connected to a pressure transducer. The device is also equipped with finger oximetry, but this feature was not used in the study. The flow signal is recorded with a sampling rate of 100 Hz and is pre-processed by linearizing the signal, filtering the noise, and zeroing the baseline. For analyzing the flow, the signal is split into epochs of 30 min as CSR occurs typically in periods of 20-30 minutes. The classifier of ApneaLink uses various features for identifying CSR. As shown in Figure 2, the main features are the cycle length (for CSR usually 30-90 s), the apneahypopnea length, the hyperpnea length, and the shape of the hyperpnea. The typical waxing and waning shape of the flow can be described in mathematical terms as a linear combination of various trigonometric functions. Furthermore, a so-called jump feature is analyzed, as displayed in Figure 3. Usually, the CSR hyperpnea starts with a moderate flow after the end of an apnea, and the hyperpnea of obstructive sleep apnea (OSA) begins typically with the so-called sudden breath. Thus, creating an envelope over a 2-s time base, the jump features of OSA and CSR are significantly different from each other. All features are combined into a pertinent set of parameters and tested







by a discriminator providing a probability of CS-like breathing between zero and one. However, if in a definite epoch the discriminator provides a high probability for CSR in patients with both CSR and OSA, it is possible that obstructive events may occur in the same epoch as well. The resulting CSR probability can be marked up on a patient record or displayed as a time-series graph. The classification of each epoch could be displayed in a bar chart (Figure 4). This provides an indication of a patient record at a glance. In this study, the results provided by ApneaLink were based on total recording time, whereas the results derived from PSG were based on total sleep time.

#### **Statistical Analysis**

Using ROC curves, we obtained the best cut-off point of ApneaLink classifier probability threshold compared to the diagnosis of sleep experts in terms of sensitivity, specificity, and positive and negative likelihood ratios.<sup>21,22</sup> The patient was considered to be positive for CSR if the probability value of CSR exceeded a specific threshold. The accuracy was the number of true positive and true negative cases divided by all cases.


#### RESULTS

The ROC curve with different probability thresholds is displayed in Figure 5. The AUC was 93.4%. Best values for sensitivity and specificity regarding the CSR probability were yielded at a threshold of 0.5. We obtained a sensitivity of 87.1% and a specificity of 94.9% as shown in Table 2. We found 2 false positive results reflecting a low CSR probability of 0.56 and 0.57 (close to the threshold of 0.5). There were 4 false nega
 Table 2—Diagnostic Accuracy for CSR Regarding Classifier

 Probability Threshold of 0.5


Sensitivity % (95% CI)	Specificity % (95% CI)	LR+	LR-	Accuracy %
87.1 (75.3-98.9)	94.9 (87.9-100)	17.0	0.14	91.4

CI denotes confidence interval, LR refers to likelihood ratio.



over a 2-s time base. (A) OSA jump feature: the maximum difference between 2 points is high compared to CSR jump feature. (B) CSR jump feature: the maximum difference between 2 points is low compared to OSA jump feature.

tive results: one included apneas of short cycle length and the other 3 contained no continuous clear CSR pattern, so that the experts relied heavily on non-flow information such as thorax and abdomen signals.




**Figure 4**—Classification of each epoch displayed in a bar chart; one epoch = 30 min. (A) A report where many epochs show a high probability of CS-like breathing. (B) A report where a single epoch shows a high probability of CS-like breathing. This epoch could be investigated more closely.

#### DISCUSSION

This study shows for the first time the utility of ApneaLink as a screening device for detecting CSR in a group of subjects with OSA, CSR, and no SDB. The ApneaLink classifier based on pattern recognition of nasal canula flow signal provided a high diagnostic accuracy for the detection of CSR when compared to the results of manually scored simultaneous PSG. The classifier of ApneaLink recognized CSR with a sensitivity of 87.1% and a specificity of 94.9%.

A limitation of ApneaLink at the current stage of development is that the classifier discriminates between CSR and non-CSR only. In theory, mixed apneas that look like central apneas of CSR are classified as CSR. Nevertheless, in this study, no false positive cases caused by incorrect classification of mixed apneas were identified.

Our study tested CSR detection by comparing ApneaLink with simultaneous PSG. Erman et al<sup>11</sup> showed that the diagnostic accuracy regarding the AHI decreased when ApneaLink in a home setting was compared with ApneaLink from an attended laboratory-based case. The reason for different diagnostic accuracy under those conditions was night-to-night variation of the occurrence of respiratory events.<sup>23,24</sup> Also, long-term loss of flow signal when the nasal cannula accidentally falls off during



sleep occurs more often in the unattended home setting than in the attended laboratory-based PSG setting. Therefore, for future purposes a study should be undertaken to evaluate the diagnostic accuracy of the CSR classifier based on results obtained from a comparison between home studies and in-laboratory studies.

To our knowledge, ApneaLink is the first commercially available screening device that is aimed to detect CSR based on pattern recognition using the flow signal. We believe that the diagnostic accuracy could be improved if the ApneaLink classifier for CSR were based on a 2-channel screening device with additional use of pulse oximetry. In a recent study El-Solh<sup>16</sup> demonstrated that a neural network based on spectral analysis of arterial oxygen saturation SpO<sub>2</sub> could detect CSR with a sensitivity of 100% and a specificity of 99%. The astonishingly high accuracy reported in this study is explained by the use of a strict definition of CSR, which was defined as >85% of SDB events being classified as central events. Based on these results, however, combining pattern recognition of flow and SpO<sub>2</sub> in a single classifier would be a promising approach for detecting CSR. The higher diagnostic accuracy of such a 2-channel device would result from a more precise classification of respiratory events if signals from one or both channels are of low quality.

In summary, this study shows that ApneaLink equipped with CSR classifier is a useful device for detecting CSR. Due to the high prevalence of CSR in patients with congestive heart failure there is a need for suitable screening devices that are easy to use and available at low cost. ApneaLink is an appropriate screening tool that may help to prioritize patients with high pre-test probability of CSR for PSG.

#### DISCLOSURE STATEMENT

This study was sponsored by a research grant from ResMed Corporation. Dr. Weinreich has participated in paid speaking engagements. Dr. Armitstead is the Manager of Engineering Research for ResMed Ltd. Dr. Töpfer has participated in paid speaking engagements. Dr. Teschler has participated in paid speaking engagements and is on the advisory board of ResMed Corporation. The other authors have indicated no financial conflicts of interest.

#### REFERENCES

- Sin DD, Fitzgerald F, Parker JD, Newton G, Floras JS, Bradley TD. Risk factors for central and obstructive sleep apnea in 450 men and women with congestive heart failure. Am J Respir Crit Care Med 1999;160:1101-6.
- 2. Cherniack NS, Longobardo G, Evangelista CJ. Causes of Cheyne-Stokes. Neurocrit Care 2005;3:271-9.
- Nopmaneejumruslers C, Kaneko Y, Hajek V, Zivanovic V, Bradley TD. Cheyne-Stokes respiration in stroke: relationship to hypocapnia and occult cardiac dysfunction. Am J Respir Crit Care Med 2005;171:1048-52.
- Oldenburg O, Lamp B, Faber L, Teschler H, Horstkotte D, Topfer V. Sleep-disordered breathing in patients with symptomatic heart failure: a contemporary study of prevalence in and characteristics of 700 patients. Eur J Heart Fail 2007;9:251-7.
- Wessendorf TE, Teschler H, Wang YM, Konietzko N, Thilmann AF. Sleep-disordered breathing among patients with first-ever stroke. J Neurol 2000;247:41-7.
- Hanly PJ, Zuberi-Khokhar NS. Increased mortality associated with Cheyne-Stokes respiration in patients with congestive heart failure. Am J Respir Crit Care Med 1996;153:272-6.
- Javaheri S, Parker TJ, Liming JD, et al. Sleep apnea in 81 ambulatory male patients with stable heart failure. Types and their prevalences, consequences, and presentations. Circulation 1998;97:2154-9.
- Javaheri S, Shukla R, Zeigler H, Wexler L. Central sleep apnea, right ventricular dysfunction, and low diastolic blood pressure are predictors of mortality in systolic heart failure. J Am Coll Cardiol 2007;49:2028-34.
- Flemons WW, Douglas NJ, Kuna ST, Rodenstein DO, Wheatley J. Access to diagnosis and treatment of patients with suspected sleep apnea. Am J Respir Crit Care Med 2004;169:668-72.
- Wang Y, Teschler T, Weinreich G, Hess S, Wessendorf TE, Teschler H. Validation of microMESAM as screening device for sleep disordered breathing. Pneumologie 2003;57:734-40.
- 11. Erman MK, Stewart D, Einhorn D, Gordon N, Casal E. Validation of the ApneaLink for the screening of sleep apnea: a novel and simple single-channel recording device. J Clin Sleep Med 2007;3:387-92.
- Bishop CM. Neural networks for pattern recognition. Oxford University Press, 1995.
- 13. Duda RO, Hart PE and Stork DG. Pattern classification. New York: Wiley, 2001.
- 14. Hastie T, Tibshirani R, Friedman J. The elements of statistical learning. Springer, 2001.
- 15. Theodoridis S, Koutroumbas K. Pattern recognition. Academic Press, 2006.
- 16. El-Solh AA, Magalang UJ, Mador MJ, et al. The utility of neu-

ral network in the diagnosis of Cheyne-Stokes respiration. J Med Eng Technol 2003 ;27:54-8.

- Weinreich G, Armitstead J, Teschler H. Pattern recognition of obstructive sleep apnoea and Cheyne-Stokes respiration. Physiol Meas 2008;29:869-78.
- Rechtschaffen A, Kales A. A Manual of standardized terminology, techniques, and scoring system for the sleep stages of human subjects. Washington, DC, U.S. Government Printing Office, 1968.
- 19. American Sleep Disorders Association. EEG arousals: scoring rules and examples. Sleep 1992;15:174-83.
- 20. American Academy of Sleep Medicine. The AASM manual for the scoring of sleep and associated events. Westchester, IL: American Academy of Sleep Medicine, 2007.

- 21. Flemons WW, Littner MR. Measuring agreement between diagnostic devices. Chest 2003;124:1535-42.
- 22. Hanley JA, McNeil BJ. The meaning and use of the area under a receiver operating characteristic (ROC) curve. Radiology 1982;143:29-36.
- Fietze I, Dingli K, Diefenbach K, et al. Night-to-night variation of the oxygen desaturation index in sleep apnoea syndrome. Eur Respir J 2004;24:987-93.
- 24. Sahlin C, Svanborg E, Stenlund H, Franklin KA. Cheyne-Stokes respiration and supine dependency. Eur Respir J 2005;25:829-33.